Puncture Pete – hyperrealistic IV training arm, light skin color





Product information "Puncture Pete – hyperrealistic IV training arm, light skin color"

Our new injection training arm **„Puncture Pete“** solves all the issues you might have with other arm. It offers a vast network of **easily replaceable veins of varying depths,** allowing for a range of difficulty levels in IV insertion practice.

The following veins can be punctured:
- dorsal arch
- cephalic vein
- basilic vein
- brachial (deep) vein
- radial vein
- median antebrachial vein
- ulnar vein
- median cubital vein

All veins that can be seen on the skin are exactly under the skin structure and can be punctured. The fluid system is capable of simulating blood flashback upon successful vein puncture, providing immediate feedback to users. The arm is **made of durable, yet realistic soft silicone** to withstand repeated use and ensure longevity in training environments. It is flexible and allows for adjustable arm positioning and rotation to facilitate varying training scenarios and user preferences. Of course it is compatible with standard IV equipment such as catheters, needles, vacuum tubes, and infusion sets.

Puncture Pete – hyperrealistic IV training arm, light skin color

Advantages:

- Skin can be changed easily by one person within some minutes. Simply roll it off and on, making the service for the arm really easy. You will love this arm as it saves you plenty of time.

- If you see a vein, you can puncture this vein! We took great lengths to ensure the surface topography corresponds to a “stickable” vessel underneath. In other words, if you can see and palpate a vein, you can stick that vein and achieve blood return.

- No more air bubbles in the veins! Air pockets lead to insufficient flashback which can leave a learner thinking they have failed, despite doing everything correctly. Our arms take less than 15 seconds to prime and purge all air bubbles from the system!

- Fast and simple vessel installation. Vessels stay put when removing/installing skin! Gone are the days of tubing getting kinked and misaligned. Retention tabs are built-in to the design of the trainer’s core to keep vessels in place during use and maintenance.

- Different skin tones – Make your training more realistic by using arms with different skin tones. The truer to life, the more real the experience - the more real the experience, the greater the transfer of knowledge.

Scope of delivery:

1x Arm, fully assembled, consisting of:

- 7710E Arm core for ‘Puncture Pete’ IV training arm
- 7710B Veins for ‘Puncture Pete’ IV training arm
- 7710C Vein distributor for ‘Puncture Pete’ IV training arm
- 7710A Skin for ‘Puncture Pete’ IV training arm, light



1x 7710D Blood bag with supply line

1x 7710F Bag for ‘Puncture Pete’ IV training arm

1x 7024 Blood-coloured liquid
1x Tripod
1x Powder

